

Contec® CritiGear™ Sterile Nitrile Gloves

USP <797> and USP <800> compliant latex-free sterile nitrile gloves

Contec CritiGear™ Sterile Nitrile Gloves have been specifically designed for pharmacy environments regulated by USP <797> and USP <800>. Clean processing and packaging during manufacturing ensures that the gloves are acceptable for ISO Class 5 applications and do not introduce harmful contaminants to products or processes.

Suitable for double gloving procedures and are appropriate for use in sterile compounding, research, and pharmaceutical manufacturing, these durable gloves have extra-long, 12 inch cuff to provide additional protection and coverage.

CritiGear Sterile Nitrile Gloves feature a textured grip and are individually pair packaged in cleanroom packaging to ensure cleanliness and cleanroom compatibility. Available in sizes X-Small to 2X-Large.

Instructions for Use: The Nitrile Examination glove is a disposable device intended for medical purpose that is worn on the examiner's hand or finger to prevent contamination between patient and examiner.

In addition, these gloves were tested for use with Chemotherapy drugs in accordance with ASTM D6978-05: Standard Practice for the Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs.



Designed for use in sterile environments and ideal for pharmacy cleanrooms. CritiGear Sterile Nitrile Gloves are also suitable for use when compounding hazardous drugs.

Features	Benefits
Tested per ASTM D6978-05	<ul style="list-style-type: none"> Protects wearer from tested chemotherapy drugs*
USP <797> compliant	<ul style="list-style-type: none"> Meets USP <797> and current draft USP <800> guidelines
Clean processed	<ul style="list-style-type: none"> Reduces harmful particulate ensuring compatibility with ISO Class 5 environments
Clean packaging	<ul style="list-style-type: none"> Pair packed in cleanroom environment to ensure glove cleanliness and cleanroom compatibility Sterile presentation packaging designed to ease proper donning of sterile gloves
Latex-free	<ul style="list-style-type: none"> Reduces potential for allergic sensitivity
Extra-long 12" length	<ul style="list-style-type: none"> Greater coverage prevents risks associated with exposed skin Protects product from potential contamination from exposed skin
Ambidextrous	<ul style="list-style-type: none"> Offers excellent fit while minimizing excessive inventory for hand-specific part numbers
Textured grip	<ul style="list-style-type: none"> Increased tactile sensitivity and dexterity
Sterility Assurance Level 10 ⁻⁶	<ul style="list-style-type: none"> Validated sterile; complete certificates of processing and sterility available upon request

Part no.	Description	Glove Size	Packaging
HCGS0002	CritiGear Sterile Nitrile Gloves, X-Small	6 - 6.5	1 pair/pack, 200 pairs/case
HCGS0012	CritiGear Sterile Nitrile Gloves, Small	7 - 7.5	1 pair/pack, 200 pairs/case
HCGS0022	CritiGear Sterile Nitrile Gloves, Medium	8 - 8.5	1 pair/pack, 200 pairs/case
HCGS0032	CritiGear Sterile Nitrile Gloves, Large	9 - 9.5	1 pair/pack, 200 pairs/case
HCGS0042	CritiGear Sterile Nitrile Gloves, X-Large	10 - 10.5	1 pair/pack, 200 pairs/case
HCGS0052	CritiGear Sterile Nitrile Gloves, 2X-Large	11	1 pair/pack, 200 pairs/case

***Notes**

a) Review chemotherapy/hazardous drug data prior to use.