

CritiGear™ Cleanroom Mask

HCCM0001

Contec® Healthcare's CritiGear™ Cleanroom Mask is a 7" ASTM F2100 Compliant Level 1 cleanroom mask, packaged in cleanroom packaging and suitable for procedures that are low risk for producing fluids, sprays, or aerosols. It is made with low-linting, lightweight thermal bond material and provides user comfort. The fluid-resistant style provides high filtration efficiency while remaining breathable. The masks are fitted with soft knit ultrasonically welded ear loops for a comfortable fit.



Medical Face Mask Material Requirements by Performance Class per ASTM F2100-11

Characteristic	Level 1
Bacterial Filtration Efficiency	% ≥ 98
Differential Pressure, mm H ₂ O/cm ²	<5.0
Sub-micron Particulate Filtration Efficiency at 0.1 micron	% ≥ 98
Resistance to Penetration by synthetic blood, minimum pressure in mm Hg for pass result	120
Flame spread	Class 1
Helmke Drum	Category 1

Test results provided by an independent third party testing laboratory.

Part No.	Description	Size	Packaging
HCCM0001	CritiGear Cleanroom Mask	One Size	50/bag; 10 bags/case